



*Innovating Health & Wellness
Within Your World*

ULTRASOUND TABLES & ACCESSORIES

- Multi-Specialty
- Echocardiography
- Vascular
- Women's Imaging
- General Ultrasound



made in the USA
with US & imported parts



ULTRASOUND TABLES

Excellent ergonomics reduce sonographer injury.

Perfect for...

- Heart & Vascular Labs
- Hospitals/Medical Centers
- Women's Health Centers
- Private Practice



Vascular with Fowler

*Inspired by the
leading ergonomic
experts in
sonography*



SOUND | ERGONOMICS

- Decrease sonographer injury
- Features minimize exam times
- 550 lb. patient weight rating

Things to consider when choosing a table:

- What types of ultrasound procedures will be performed currently and in the future?
- What electric movements are a must for a new table?
- Will the table be moved around the room or mainly positioned in one spot?
- What accessories will be required?
- Are there any ergonomic issues to be addressed with the purchase of a new table?
- What is the budget?



**OAKWORKS
MEDICAL**



See available models
on pages 11-13

Start with a General Ultrasound Table, then simply add the features you want based on the type of ultrasound specialty you perform.

- Create the perfect table for what you do
- Don't pay for features you don't want
- Up to 12 different options to choose from

*Build-A-Bed™ is Trademarked
by Sound Ergonomics.*

Benefits...

For The Sonographer

- Build-A-Bed concept enables you to customize the table with all the features you need
- Design supports good body mechanics, minimizing strain and injury.
- Unobstructive wheel base design allows mobility around the table.

For The Patient

- Ultimate in comfort with our 3" padding
- Easy on & off, wheelchair accessible
- 550 lb. bariatric rated lift & load capacity

For The Facility

- Customizing your table is easy & cost effective, only pay for the features you want
- Low maintenance, all electric design
- 3 years parts, 2 years labor warranty

Choose a Specialty Table

These models were designed to meet the needs of almost every kind of facility. You can add other options to create the perfect table for your specific requirements.

- 5 preconfigured models to choose from.
- Features grouped to maximize your efficiency.
- Add options to create your perfect table.



Multi-Specialty

PRODUCTS AT A GLANCE

Multi-Specialty Table

#67375-T* 27" W #61803-T* 30" W



FEATURES

Top: 27" (69cm) or 30" (76cm) wide x 72" (183cm) long
Electric Height Range: 22" - 38" (56-97cm)
Electric Trendelenburg/Reverse Tilt: $\pm 15^\circ$
Electric Fowler Backrest 0-75°
Manual 4 Position Leg Rest 0-90°
Hand Control & Foot Control
Weight Capacity: 550 lb (250kg)
Padding: 3" (8cm) Comfort Foam™
4" (10cm) Locking Casters
Adjustable Access/Back Support - Left 0-85°
Adjustable Access/Back Support - Right 0-85°
Central Base Locking System
Paper Roll Holder
Stirrups
Warranty: 3 years parts, 2 years labor
FDA listed, cETLus listed, CB Certificate

OPTIONS *must be ordered with table*

	#
Accessory Mounting Plate	62123
Foot Control (Height, Tilt, Fowler, Cardio Panel)	65679
Non-USA electric system (contact us for details)	73134
Procedure Tray (30" W Only)	73127
Side Rails (set of 2)	61781
Sliding Cardio Panel - Electric (30" W, Left Side)	73131
Sliding Cardio Panel - Manual (30" W, Left Side)	73130

Echocardiography Table

#65258-T*



FEATURES

Top: 30" (76cm) wide x 72" (183cm) long
Electric Height Range: 22" - 38" (56-97cm)
Electric Trendelenburg/Reverse Tilt: $\pm 15^\circ$
Electric Fowler Backrest 0-75°
Hand Control
Weight Capacity: 550 lb (250kg)
Padding: 3" (8cm) Comfort Foam™
4" (10cm) Locking Casters
Adjustable Access/Back Support - Right 0-85°
Manual Sliding Cardio Panel - Left Side
Paper Roll Holder
Warranty: 3 years parts, 2 years labor
FDA listed, cETLus listed, CB Certificate

OPTIONS *must be ordered with table*

	#
Accessory Mounting Plate	62123
Central Base Locking System (<i>can not be used with Patient Lift Compatible Base</i>)	61777
Foot Control (Height, Tilt, Fowler)	61901
Foot Control (Height, Tilt, Fowler, Cardio Panel)	65679
Non-USA electric system (contact us for details)	73134
Patient Lift Compatible Base (<i>can not be used with Central Base Locking System</i>)	65893
Side Rails (set of 2)	61781
Sliding Cardio Panel - Electric Upgrade	73131

Vascular Table with Fowler

#66783-T* 27" W #64219-T* 30" W



FEATURES

Top: 27" (69cm) or 30" (76cm) wide x 72" (183cm) long
Electric Height Range: 23" - 40" (58-102cm)
Electric Trendelenburg 15°
Reverse Trendelenburg: 40° + Electric Fowler = 70°
Electric Fowler Backrest
Hand Control
Weight Capacity: 550 lb (250kg)
Padding: 3" (8cm) Comfort Foam™
4" (10cm) Locking Casters
Integrated Foot Board (split, drop away)
Paper Roll Holder
Patient Safety Strap
Side Rails (set of 2)
Warranty: 3 years parts, 2 years labor
FDA listed, cETLus listed, CB Certificate

OPTIONS *must be ordered with table*

	#
Accessory Mounting Plate	62123
Adjustable Access/Back Support - Left 0-85°	73133
Adjustable Access/Back Support - Right 0-85°	73132
Central Base Locking System (<i>can not be used with Patient Lift Compatible Base</i>)	61777
Foot Control (Height, Tilt, Fowler, Vascular)	65679
Non-USA electric system (contact us for details)	73134
Patient Lift Compatible Base (<i>can not be used with Central Base Locking System</i>)	65893
Sliding Cardio Panel - Manual (30" W, Left Side)	73130

Vascular Table without Fowler

#70217-T* 27" W #64928-T* 30" W



FEATURES

Top: 27" (69cm) or 30" (76cm) wide x 72" (183cm) long
Electric Height Range: 23" - 40" (58-102cm)
Electric Trendelenburg 15°
Reverse Trendelenburg: 40°
Hand Control
Weight Capacity: 550 lb (250kg)
Padding: 3" (8cm) Comfort Foam™
4" (10cm) Locking Casters
Integrated Foot Board (split, drop away)
Paper Roll Holder
Patient Safety Strap
Side Rails (set of 2)
Warranty: 3 years parts, 2 years labor
FDA listed, cETLus listed, CB Certificate

OPTIONS *must be ordered with table*

	#
Accessory Mounting Plate	62123
Adjustable Access/Back Support - Left 0-85°	73133
Adjustable Access/Back Support - Right 0-85°	73132
Central Base Locking System (<i>can not be used with Patient Lift Compatible Base</i>)	61777
Foot Control (Height, Tilt, Vascular)	61901
Non-USA electric system (contact us for details)	73134
Patient Lift Compatible Base (<i>can not be used with Central Base Locking System</i>)	65893
Sliding Cardio Panel - Manual (30" W, Left Side)	73130



OAKWORKS
MEDICAL

See options and accessories on page 14-15.

*Fabric color code must be specified. Fabric options on back cover.

PRODUCTS AT A GLANCE

Women's Imaging Table

#66127-T* 27" W #61848-T* 30" W



FEATURES

Top: 27" (69cm) or 30" (76cm) wide x 72" (183cm) long
Electric Height Range: 22" - 38" (56-97cm)
Electric Trendelenburg/Reverse Tilt: $\pm 15^\circ$
Electric Fowler Backrest 0-75°
Manual 4 Position Leg Rest 0-90°
Hand Control
Weight Capacity: 550 lb (250kg)
Padding: 3" (8cm) Comfort Foam™
4" (10cm) Locking Casters
Adjustable Access/Back Support - Right 0-85°
Paper Roll Holder
Stirrups
Warranty: 3 years parts, 2 years labor
FDA listed, cETLus listed, CB Certificate

OPTIONS *must be ordered with table*

	#
Accessory Mounting Plate	62123
Adjustable Access/Back Support - Left 0-85°	73133
Central Base Locking System (can not be used with Patient Lift Compatible Base)	61777
Foot Control (Height, Tilt, Fowler)	61901
Foot Control (Height, Tilt, Fowler, Cardio Panel)	65679
Non-USA electric system (contact us for details)	73134
Patient Lift Compatible Base (can not be used with Central Base Locking System)	65893
Procedure Tray (30" W Only)	73127
Side Rails (set of 2)	61781
Sliding Cardio Panel - Electric (30" W, Left Side)	74080
Sliding Cardio Panel - Manual (30" W, Left Side)	73130

**BUILD
-A-
BED™**

General Rectangular Top Table

#68666-T* 27" W #61989-T* 30" W



FEATURES

Top: 27" (69cm) or 30" (76cm) wide x 72" (183cm) long
Electric Height Range: 22" - 38" (56-97cm)
Electric Trendelenburg/Reverse Tilt: $\pm 15^\circ$
Hand Control
Weight Capacity: 550 lb (250kg)
Padding: 3" (8cm) Comfort Foam™
4" (10cm) Locking Casters
Paper Roll Holder
Warranty: 3 years parts, 2 years labor
FDA listed, cETLus listed, CB Certificate

OPTIONS *must be ordered with table*

	#
Accessory Mounting Plate	62123
Adjustable Access/Back Support - Left 0-85°	73133
Adjustable Access/Back Support - Right 0-85°	73132
Central Base Locking System (can not be used with Patient Lift Compatible Base)	61777
Foot Control (Height, Tilt)	61901
Non-USA electric system (contact us for details)	73134
Patient Lift Compatible Base (can not be used with Central Base Locking System)	65893
Side Rails (set of 2)	61781
Sliding Cardio Panel - Manual (30" W, Left Side)	73130

**BUILD
-A-
BED™**

General Table with Fowler

#69222-T* 27" W #61988-T* 30" W



FEATURES

Top: 27" (69cm) or 30" (76cm) wide x 72" (183cm) long
Electric Height Range: 22" - 38" (56-97cm)
Electric Trendelenburg/Reverse Tilt: $\pm 15^\circ$
Electric Fowler Backrest 0-75°
Hand Control
Weight Capacity: 550 lb (250kg)
Padding: 3" (8cm) Comfort Foam™
4" (10cm) Locking Casters
Paper Roll Holder
Warranty: 3 years parts, 2 years labor
FDA listed, cETLus listed, CB Certificate

OPTIONS *must be ordered with table*

	#
Accessory Mounting Plate	62123
Adjustable Access/Back Support - Left 0-85°	73133
Adjustable Access/Back Support - Right 0-85°	73132
Central Base Locking System (can not be used with Patient Lift Compatible Base)	61777
Foot Control (Height, Tilt, Fowler)	61901
Foot Control (Height, Tilt, Fowler, Cardio Panel)	65679
Non-USA electric system (contact us for details)	73134
Patient Lift Compatible Base (can not be used with Central Base Locking System)	65893
Side Rails (set of 2)	61781
Sliding Cardio Panel - Electric (30" W, Left Side)	74080
Sliding Cardio Panel - Manual (30" W, Left Side)	73130

**BUILD
-A-
BED™**

General 3-Section Top Table

#62963-T* 27" W #63028-T* 30" W



FEATURES

Top: 27" (69cm) or 30" (76cm) wide x 72" (183cm) long
Electric Height Range: 22" - 38" (56-97cm)
Electric Trendelenburg/Reverse Tilt: $\pm 15^\circ$
Electric Fowler Backrest 0-75°
Manual 4 Position Leg Rest 0-90°
Hand Control
Weight Capacity: 550 lb (250kg)
Padding: 3" (8cm) Comfort Foam™
4" (10cm) Locking Casters
Paper Roll Holder
Warranty: 3 years parts, 2 years labor
FDA listed, cETLus listed, CB Certificate

OPTIONS *must be ordered with table*

	#
Accessory Mounting Plate	62123
Adjustable Access/Back Support - Left 0-85°	73133
Adjustable Access/Back Support - Right 0-85°	73132
Central Base Locking System (can not be used with Patient Lift Compatible Base)	61777
Foot Control (Height, Tilt, Fowler)	61901
Foot Control (Height, Tilt, Fowler, Cardio Panel)	65679
Non-USA electric system (contact us for details)	73134
Patient Lift Compatible Base (can not be used with Central Base Locking System)	65893
Procedure Tray (30" W Only)	73127
Side Rails (set of 2)	61781
Sliding Cardio Panel - Electric (30" W, Left Side)	74080
Sliding Cardio Panel - Manual (30" W, Left Side)	73130
Stirrups	71001



**OAKWORKS
MEDICAL**

See options and accessories on page 14-15.

*Fabric color code must be specified. Fabric options on back cover.

MULTI-SPECIALTY TABLE

- 550 lb (250kg) bariatric rated lift & load capacity
- Twin tower design for added strength & rigidity
- Dynamic 22" to 38" (56-97cm) adjustable height range for wheelchair transfer & ergonomic working heights



Easily converts to fully upright chair position in seconds!

Recommended for:

- All abdominal & superficial structures exams (thyroid, male pelvis)
- Musculoskeletal exams
- OB/GYN exams
- Breast exams
- Vascular exams
- Biopsies



4 Position Leg Section

MULTI-SPECIALTY TABLE

#67375-T* 27" W

#61803-T* 30" W

FEATURES	
3 Section Top	27" (69cm) or 30" (76cm) wide x 72" (183cm) long
Backrest	34" (86cm) long
Seat	24" (61cm) long
Leg Rest	13" (33cm) long
Motion 1	Electric Height Range 22" - 38" (56-97cm)
Motion 2	Electric Trendelenburg/Reverse Tilt $\pm 15^\circ$
Motion 3	Electric Fowler Backrest 0-75°
Motion 4	Manual 4 position Leg Rest 0-90°
Controls	Multi-function hand & foot control
Weight Capacity	550 lb (250kg)
Product Weight	359 lb (163kg)
Padding	3" (8cm) Comfort Foam™
Included Features	4" (10cm) locking casters, Adjustable access back support (left & right) 0-85°, Central base locking system, Paper roll holder, Stirrups (Retractable foot rests)
Warranty	3 years parts, 2 years labor
Safety Listings	FDA listed, cETLus listed, CB Certificate

OPTIONS <i>must be ordered with table. See page 14-15 for details.</i>		#
Accessory Mounting Plate	Fowler end	62123
Foot Control	Height, Tilt, Fowler, Cardio Panel	65679
Non-USA electric system	please contact us for order details	73134
Procedure Tray 30" (76 cm) Width Only	12.25" (6cm) x 10.25" (26cm) and holds 2 quarts (1.9 liters)	73127
Side Rails (set of 2)	23.5" (60cm) long x 6" (15cm) tall (from surface of table)	61781
Sliding Cardio Panel - Electric	30" (76cm) width table, left side	73131
Sliding Cardio Panel - Manual	30" (76cm) width table, left side	73130
ACCESSORIES <i>order any time. See accessories on page 15.</i>		#
Detachable Arm Board	Use with T-rail adapter 5" (13cm) wide x 20.5" (52cm) long	62192
QuickLock™ Head Rest (black) & AeroCel Pad	Comfortable, adjustable prone or supine support (requires mounting plate)	PKG6914-T*
Table Extender	24" (61cm) wide x 12" (30cm) long (requires mounting plate)	63487
T-rail Adapter 8" (20cm) long	Attaches to the backrest for use with USA format T-rail attachable accessories.	61902

*Fabric color code must be specified. Fabric options on back cover.

ECHOCARDIOGRAPHY TABLE

- Sliding cardio panel allows you to adjust the opening to the size you need
- Electric or manual sliding cardiac panel opens as wide as 15" (38cm). Great for left or right handed scanning
- 22" to 38" (56-97cm) height range supports ergonomic working heights for both echocardiograms & stress echo
- Designed specifically for vascular & cardiac exams
- Dual tower design for added strength & rigidity allowing a 550 lb (250kg) pound lift/load capacity



Adjustable Access/Back Support



Sliding Cardio Panel opens as wide as 15" (38cm)

Recommended for:

- Cardiac exams
- Vascular exams
- All abdominal & superficial structures exams
- Musculoskeletal exams
- Biopsies

ECHOCARDIOGRAPHY TABLE

#65258-T*

FEATURES	
2 Section Top	30" (76cm) wide x 72" (183cm) long
Backrest	34" (86cm) long
Seat	38" (97cm) long
Motion 1	Electric Height Range 22" - 38" (56-97cm)
Motion 2	Electric Trendelenburg/Reverse Tilt $\pm 15^\circ$
Motion 3	Electric Fowler Backrest 0-75°
Motion 4	Manual Sliding Cardio Panel opens 15" (38cm)
Controls	Multi-function hand control
Weight Capacity	550 lb (250kg)
Product Weight	359 lb (163kg)
Padding	3" (8cm) Comfort Foam™
Included Features	4" (10cm) locking casters, Adjustable access back support (right) 0-85°, Manual sliding cardio panel (left) opens 15" (38cm), Paper roll holder
Warranty	3 years parts, 2 years labor
Safety Listings	FDA listed, cETLus listed, CB Certificate

OPTIONS must be ordered with table. See page 14-15 for details.		#
Accessory Mounting Plate	Fowler end	62123
Central Base Locking System	Not available with Patient Lift Compatible Base	61777
Foot Control	Height, Trendelenburg, Backrest	61901
Foot Control	Height, Trendelenburg, Backrest, Cardio Panel	65679
Non-USA electric system	please contact us for order details	73134
Patient Lift Compatible Base	Not available with Central Base Locking System	65893
Side Rails (set of 2)	23.5" (60cm) long x 6" (15cm) tall (from surface of table)	61781
Sliding Cardio Panel - Electric	Upgrade your manual panel opens 15" (38cm)	73131
ACCESSORIES order any time. See accessories on page 15.		#
Detachable Arm Board	Use with T-rail adapter 5" (13cm) wide x 20.5" (52cm) long	62192
T-rail Adapter 8" (20cm) long	Attaches to the backrest for use with USA format T-rail attachable accessories.	61902
Table Extender	24" (61cm) wide x 12" (30cm) long (requires mounting plate)	63487
QuickLock™ Head Rest (black) & AeroCel Pad	Comfortable, adjustable prone or supine support (requires mounting plate)	PKG6914-T*

*Fabric color code must be specified. Fabric options on back cover.

VASCULAR TABLE WITH FOWLER

- Dual towers & electric fowler makes 70° reverse Trendelenburg a reality†
- Safely position clients for Venous Reflux studies.
- Split foot board design allows non-weight bearing exams and fold down out of the way for carotid artery exams. (see page 9)



Easy Hand Control Adjustments
Up To 70° Reverse Trendelenburg



Recommended for:

- Venous reflux studies
- All other vascular exams
- Abdominal & superficial structures exams
- Musculoskeletal exams
- Biopsies

“Finally we have a table designed for vascular ultrasound!”

Patricia A Poe, BA RVT FSVU

Director Noninvasive Vascular Laboratory Services,
UC Health

VASCULAR TABLE WITH FOWLER

#66783-T* 27" W

#64219-T* 30" W

FEATURES	
2 Section Top	27" (69cm) or 30" (76cm) wide x 72" (183cm) long
Backrest	34" (86cm) long
Seat	38" (97cm) long
Motion 1	Electric Height Range 23" - 40" (58-102cm)
Motion 2	Electric Trendelenburg/Reverse Tilt $\pm 15^\circ$
Motion 3	Electric Additional Reverse Trendelenburg 0-25° (40° total)
Motion 4†	Electric 0-75° Backrest (at 40° max tilt, 70° functional)
Controls	Multi-function hand control
Weight Capacity	550 lb (250kg)
Product Weight	452 lb (205kg)
Padding	3" (8cm) Comfort Foam™
Included Features	4" (10cm) locking casters, Integrated foot board (split, drop away), Paper roll holder, Patient safety strap, Side rails (2) 23.5" (60cm) long x 6" (15cm) tall (from surface of table)
Warranty	3 years parts, 2 years labor
Safety Listings	FDA listed, cETLus listed, CB Certificate

OPTIONS must be ordered with table. See page 14-15 for details.		#
Accessory Mounting Plate	Fowler end	62123
Adjustable Access/Back Support	Left side 0-85°	73133
Adjustable Access/Back Support	Right side 0-85°	73132
Central Base Locking System	Not available with Patient Lift Compatible Base	61777
Foot Control	Height, Trendelenburg, Backrest, Vascular	65679
Non-USA electric system	Please contact us for order details	73134
Patient Lift Compatible Base	Not available with Central Base Locking System	65893
Sliding Cardio Panel - Manual	30" (76cm) width table, left side opens 15" (38cm)	73130
ACCESSORIES order any time. See accessories on page 15.		#
Detachable Arm Board	Use with T-rail adapter 5" (13cm) wide x 20.5" (52cm) long	62192
QuickLock™ Head Rest (black) & AeroCel Pad	Comfortable, adjustable prone or supine support (requires mounting plate)	PKG6914-T*
Table Extender	24" (61cm) wide x 12" (30cm) long (requires mounting plate)	63487
T-rail Adapter 8" (20cm) long	Attaches to the backrest for use with USA format T-rail attachable accessories.	61902

*Fabric color code must be specified. Fabric options on back cover.

VASCULAR TABLE WITHOUT FOWLER

- Dynamic 23" to 40" (58-102cm) adjustable height range for wheelchair transfer & ergonomic working heights
- Angles up to 40°
- Split foot board design allows non-weight bearing exams and fold down out of the way for carotid artery exams.
- Plenty of room for your legs or equipment under the table



Integrated Foot Boards Allow Staggered Positions, Drop Out Of The Way & Can Be Removed



Recommended for:

- All vascular **except** reflux studies
- Abdominal & superficial structures exams
- Musculoskeletal exams
- Biopsies

VASCULAR TABLE #70217-T* 27" W #64928-T* 30" W

FEATURES	
Rectangular Top	27" (69cm) or 30" (76cm) wide x 72" (183cm) long
Motion 1	Electric Height Range 23" - 40" (58-102cm)
Motion 2	Electric Trendelenburg/Reverse Tilt $\pm 15^\circ$
Motion 3	Electric Additional Reverse Trendelenburg 0-25° (40° total)
Controls	Multi-function hand control
Weight Capacity	550 lb (250kg)
Product Weight	452 lb (205kg)
Padding	3" (8cm) Comfort Foam™
Included Features	4" (10cm) locking casters, Integrated foot board (split, drop away), Paper roll holder, Patient safety strap, Side rails (2) 23.5" (60cm) long x 6" (15cm) tall (from surface of table)
Warranty	3 years parts, 2 years labor
Safety Listings	FDA listed, cETLus listed, CB Certificate

OPTIONS <i>must be ordered with table. See page 14-15 for details.</i>		#
Accessory Mounting Plate	Fowler end	62123
Adjustable Access/Back Support	Left side 0-85°	73133
Adjustable Access/Back Support	Right side 0-85°	73132
Central Base Locking System	Not available with Patient Lift Compatible Base	61777
Foot Control	Height, Trendelenburg, Vascular	61901
Non-USA electric system	Please contact us for order details	73134
Patient Lift Compatible Base	Not available with Central Base Locking System	65893
Sliding Cardio Panel - Manual	30" (76cm) width table, left side opens 15" (38cm)	73130
ACCESSORIES <i>order any time. See accessories on page 15.</i>		#
Detachable Arm Board	Use with T-rail adapter 5" (13cm) wide x 20.5" (52cm) long	62192
QuickLock™ Head Rest (black) & AeroCel Pad	Comfortable, adjustable prone or supine support (requires mounting plate)	PKG6914-T*
Table Extender	24" (61cm) wide x 12" (30cm) long (requires mounting plate)	63487
T-rail Adapter 8" (20cm) long	Attaches to the backrest for use with USA format T-rail attachable accessories.	61902

*Fabric color code must be specified. Fabric options on back cover.

WOMEN'S IMAGING TABLE

- Full featured 3 section table top design for all Women's Imaging ultrasound modalities
- 550 lb (250kg) bariatric rated lift & load capacity
- Twin tower design for added strength and rigidity
- Sleek base design promotes unobstructed compatibility with ultrasound equipment ensuring better ergonomic conditions for the sonographer



Adjustable Access Back Support



Stirrups

Recommended for:

- OB/GYN exams
- All abdominal & superficial structures exams
- Musculoskeletal exams
- Vascular exams
- Biopsies

WOMEN'S IMAGING TABLE

#66127-T* 27" W #61848-T* 30" W

FEATURES	
3 Section Top	27" (69cm) or 30" (76cm) wide x 72" (183cm) long
Backrest	34" (86cm) long
Seat	24" (61cm) long
Leg Rest	13" (33cm) long
Motion 1	Electric Height Range 22" - 38" (56-97cm)
Motion 2	Electric Trendelenburg/Reverse Tilt $\pm 15^\circ$
Motion 3	Electric Fowler Backrest 0-75°
Motion 4	Manual 4 position Leg Rest 0-90°
Controls	Multi-function hand control
Weight Capacity	550 lb (250kg)
Product Weight	359 lb (163kg)
Padding	3" (8cm) Comfort Foam™
Included Features	4" (10cm) locking casters, Adjustable access back support (right) 0-85°, Paper roll holder, Stirrups (Retractable foot rests)
Warranty	3 years parts, 2 years labor
Safety Listings	FDA listed, cETLus listed, CB Certificate

OPTIONS must be ordered with table. See page 14-15 for details.		#
Accessory Mounting Plate	Fowler end	62123
Adjustable Access/Back Support	Left side 0-85°	73133
Foot Control	Height, Trendelenburg, Fowler	61901
Foot Control	Height, Trendelenburg, Fowler, Cardio Panel	65679
Central Base Locking System	Not available with Patient Lift Compatible Base	61777
Non-USA electric system	please contact us for order details	73134
Patient Lift Compatible Base	Not available with Central Base Locking System	65893
Procedure Tray 30" (76 cm) Width Only	12.25" (6cm) x 10.25" (26cm) and holds 2 quarts (1.9 liters)	73127
Side Rails (set of 2)	23.5" (60cm) long x 6" (15cm) tall (from surface of table)	61781
Sliding Cardio Panel - Electric	30" (76cm) width table, left side	74080
Sliding Cardio Panel - Manual	30" (76cm) width table, left side	73130
ACCESSORIES order any time. See accessories on page 15.		#
Detachable Arm Board	Use with T-rail adapter 5" (13cm) wide x 20.5" (52cm) long	62192
QuickLock™ Head Rest (black) & AeroCel Pad	Comfortable, adjustable prone or supine support (requires mounting plate)	PKG6914-T*
Table Extender	24" (61cm) wide x 12" (30cm) long (requires mounting plate)	63487
T-rail Adapter 8" (20cm) long	Attaches to the backrest for use with USA format T-rail attachable accessories.	61902

*Fabric color code must be specified. Fabric options on back cover.

GENERAL RECTANGULAR TOP TABLE

- Flat top table with many optional features to customize your needs
- 550 lb (250kg) bariatric rated lift and load capacity



Trendelenburg/Reverse



22" - 38" Height Range

Recommended for:

- Abdominal & superficial structures exams
- Musculoskeletal exams
- Vascular exams
- Biopsies

GENERAL RECTANGULAR TOP

#68666-T* 27" W #61989-T* 30" W

FEATURES	
Rectangular Top	27" (69cm) or 30" (76cm) wide x 72" (183cm) long
Motion 1	Electric Height Range 22" - 38" (56-97cm)
Motion 2	Electric Trendelenburg/Reverse Tilt $\pm 15^\circ$
Controls	Multi-function hand control
Weight Capacity	550 lb (250kg)
Product Weight	359 lb (163kg)
Padding	3" (8cm) Comfort Foam™
Fabric	Premium, PVC free TerraTouch™
Included Features	4" (10cm) locking casters, Paper roll holder
Warranty	3 years parts, 2 years labor
Safety Listings	FDA listed, cETLus listed, CB Certificate

OPTIONS <i>must be ordered with table. See page 14-15 for details.</i>		#
Accessory Mounting Plate	Either end	62123
Adjustable Access Back Support	Left side 0-85°	73133
Adjustable Access Back Support	Right side 0-85°	73132
Central Base Locking System	Not available with Patient Lift Compatible Base	61777
Foot Control	Height, Trendelenburg	61901
Non-USA electric system	Please contact us for order details	73134
Patient Lift Compatible Base	Not available with Central Base Locking System	65893
Side Rails (set of 2)	23.5" (60cm) long x 6" (15cm) tall (from surface of table)	61781
Sliding Cardio Panel - Manual	30" (76cm) width table, left side opens 15" (38cm)	73130
ACCESSORIES <i>order any time. See accessories on page 15.</i>		#
Detachable Arm Board	Use with T-rail adapter 5" (13cm) wide x 20.5" (52cm) long	62192
QuickLock™ Head Rest (black) & AeroCel Pad	Comfortable, adjustable prone or supine support (requires mounting plate)	PKG6914-T*
Table Extender	24" (61cm) wide x 12" (30cm) long (requires mounting plate)	63487
T-rail Adapter 8" (20cm) long	Attaches to the backrest for use with USA format T-rail attachable accessories.	61902

*Fabric color code must be specified. Fabric options on back cover.

GENERAL TABLE WITH FOWLER

- 2 section top table with many optional features to customize your needs
- Dynamic 22" to 38" (56-97cm) adjustable height range for wheelchair transfer & ergonomic working heights
- 550 lb (250kg) bariatric rated lift & load capacity
- Twin tower design for added strength & rigidity



Paper Roll Holder folds out of the way



Hand Control

GENERAL TABLE WITH FOWLER

#69222-T* 27" W #61988-T* 30" W

FEATURES	
2 Section Top	27" (69cm) or 30" (76cm) wide x 72" (183cm) long
Backrest	34" (86cm) long
Seat	38" (97cm) long
Motion 1	Electric Height Range 22" - 38" (56-97cm)
Motion 2	Electric Trendelenburg/Reverse Tilt $\pm 15^\circ$
Motion 3	Electric Fowler Backrest 0-75°
Controls	Multi-function hand control
Weight Capacity	550 lb (250kg)
Product Weight	359 lb (163kg)
Padding	3" (8cm) Comfort Foam™
Fabric	Premium, PVC free TerraTouch™
Included Features	4" (10cm) locking casters, Paper roll holder
Warranty	3 years parts, 2 years labor
Safety Listings	FDA listed, cETLus listed, CB Certificate

OPTIONS must be ordered with table. See page 14-15 for details.		#
Accessory Mounting Plate	Fowler end	62123
Adjustable Access Back Support	Left side 0-85°	73133
Adjustable Access Back Support	Right side 0-85°	73132
Central Base Locking System	Not available with Patient Lift Compatible Base	61777
Foot Control	Height, Trendelenburg, Backrest	61901
Foot Control	Height, Trendelenburg, Backrest, Cardio Panel	65679
Non-USA electric system	please contact us for order details	73134
Patient Lift Compatible Base	Not available with Central Base Locking System	65893
Side Rails (set of 2)	23.5" (60cm) long x 6" (15cm) tall (from surface of table)	61781
Sliding Cardio Panel - Electric	30" (76cm) width table, left side opens 15" (38cm)	73131
Sliding Cardio Panel - Manual	30" (76cm) width table, left side opens 15" (38cm)	73130
ACCESSORIES order any time. See accessories on page 15.		#
Detachable Arm Board	Use with T-rail adapter 5" (13cm) wide x 20.5" (52cm) long	62192
QuickLock™ Head Rest (black) & AeroCel Pad	Comfortable, adjustable prone or supine support (requires mounting plate)	PKG6914-T*
Table Extender	24" (61cm) wide x 12" (30cm) long (requires mounting plate)	63487
T-rail Adapter 8" (20cm) long	Attaches to the backrest for use with USA format T-rail attachable accessories.	61902

*Fabric color code must be specified. Fabric options on back cover.

GENERAL 3 SECTION TOP TABLE

- Three section top table with 12 optional features to customize your table
- Sleek base design promotes unobstructed compatibility with ultrasound equipment ensuring better ergonomic conditions for the sonographer



Easily Converts To Fully Upright Chair Position In Seconds!



4 Position Leg Section



Recommended for:

- All abdominal & superficial structures exams
- OB/GYN exams
(limited use with endovaginal exams)
- Musculoskeletal exams
- Vascular exams
- Biopsies

GENERAL 3 SECTION TOP TABLE

#62963-T* 27" W #63028-T* 30" W

FEATURES	
3 Section Top	27" (69cm) or 30" (76cm) wide x 72" (183cm) long
Backrest	34" (86cm) long
Seat	24" (61cm) long
Leg Rest	13" (33cm) long
Motion 1	Electric Height Range 22" - 38" (56-97cm)
Motion 2	Electric Trendelenburg/Reverse Tilt $\pm 15^\circ$
Motion 3	Electric Fowler Backrest 0-75°
Motion 4	Manual 4 position Leg Rest 0-90°
Controls	Multi-function hand control
Weight Capacity	550 lb (250kg)
Product Weight	359 lb (163kg)
Padding	3" (8cm) Comfort Foam™
Included Features	4" (10cm) locking casters, Paper roll holder
Warranty	3 years parts, 2 years labor
Safety Listings	FDA listed, cETLus listed, CB Certificate

OPTIONS must be ordered with table. See page 14-15 for details.		#
Accessory Mounting Plate	Fowler end	62123
Adjustable Access/Back Support	Left side 0-85°	73133
Adjustable Access/Back Support	Right side 0-85°	73132
Central Base Locking System	Not available with Patient Lift Compatible Base	61777
Foot Control	Height, Trendelenburg, Fowler	61901
Foot Control	Height, Trendelenburg, Fowler, Cardio Panel	65679
Non-USA electric system	please contact us for order details	73134
Patient Lift Compatible Base	Not available with Central Base Locking System	65893
Procedure Tray	12.25" (6cm) x 10.25" (26cm) and holds 2 quarts (1.9 liters)	73127
Side Rails (set of 2)	23.5" (60cm) long x 6" (15cm) tall (from surface of table)	61781
Sliding Cardio Panel - Electric	30" (76cm) width table, left side	74080
Sliding Cardio Panel - Manual	30" (76cm) width table, left side	73130
Stirrups (Retractable Foot Rests)	Pull out from under table, fully adjustable	71001
ACCESSORIES order any time. See accessories on page 15.		#
Detachable Arm Board	Use with T-Rail adapter 5" (13cm) wide x 20.5" (52cm) long	62192
QuickLock™ Head Rest (black) & AeroCel Pad	Comfortable, adjustable prone or supine support (requires mounting plate)	PKG6914-T*
Table Extender	24" (61cm) wide x 12" (30cm) long (requires mounting plate)	63487
T-Rail Adapter 8" (20cm) long	Attaches to the backrest for use with USA format T-Rail attachable accessories.	61902

*Fabric color code must be specified. Fabric options on back cover.

ULTRASOUND FEATURES & OPTIONS

ADJUSTABLE ACCESS BACK SUPPORT

Left: #73133 Right: #73132

Adjustable with 5 positions 0-85°, allowing easy access for sonographer to the far side & back support for patients.
10" (25cm) wide x 15.5" (39cm) long

Standard Feature: Right Side - Multi-Specialty, Echocardiography, & Women's Imaging
Left Side - Multi-Specialty

Option: Right Side - Vascular Table with Fowler, Vascular Table without Fowler, General 3-Section, General with Fowler, & General Rectangular Top

Left Side - Vascular Table with Fowler, Vascular Table without Fowler, Women's Imaging, General 3-Section, General with Fowler, & General Rectangular Top



SLIDING CARDIO PANEL *patent pending*

Electric: #73131 Manual: #73130

The sliding cardio panel allows you to adjust the opening to the size you need, up to 15" (38cm). 10" (25cm) wide x 28" (71cm) long. Great for left or right handed scanning. Can also be used as an armrest. Available in manual or electric adjustment options. Available on left side of 30" (76cm) width tables only.

Standard Feature: Manual - Echocardiography

Option: Manual - Multi-Specialty, Vascular Table with Fowler, Vascular Table without Fowler, Women's Imaging, General 3-Section, General with Fowler, & General Rectangular Top

Electric - Multi-Specialty, Echocardiography, Women's Imaging, General 3-Section, & General with Fowler



PATIENT LIFT COMPATIBLE BASE #65893

This option allows the use of patient lifts & easy wheel chair transfers to minimize or eliminate caregiver injuries associated with patient transfers. The high clearance base permits all patient lifts to slide easily under the center of the table for stable transfers while not increasing the low height range of the tables that enable easy wheel chair transfers. Must be ordered with table.

- Includes larger, 5" (13cm) diameter Heavy Duty locking casters
- Provides 6.5" (17cm) clearance for patient lift
- Central Base Locking System not available with this option

Option: Echocardiography, Vascular Table with Fowler, Vascular Table without Fowler, Women's Imaging, General 3-Section, General with Fowler, & General Rectangular Top



FULL FUNCTION CONTROLS

#enquire for part number

You may choose to have two controls: two foot or one hand & one foot.

Option: All Exam Tables & Procedure Chairs



Multi-Function Foot Control



Multi-Function Hand Control

PROCEDURE TRAY #73127

Retractable, removable stainless steel. 12.25" (6cm) x 10.25" (26cm) & holds 2 quarts (1.9 liters). Available on 30" (76cm) wide tables only.

Option: Multi-Specialty, Women's Imaging, & General 3-Section



ACCESSORY MOUNTING PLATE

#62123

Add a face rest or table extender to your table.

Option: all Ultrasound tables.



*Fabric color code must be specified. Fabric options on back cover.

ULTRASOUND FEATURES & OPTIONS

STIRRUPS #71001

Adjust to a wide range of angles & lengths to accommodate any patient or procedure.

Standard Feature: Multi-Specialty & Women's Imaging

Option: General 3-Section



SIDE RAILS #61781

Increase stability & security for your patient. Fold down out of the way when not in use. 23.5" (60cm) long x 6" (15cm) tall (from surface of table)

Standard Feature: Vascular Tables

Option: Multi-Specialty, Echocardiography, Women's Imaging, General 3-Section, General with Fowler, General Rectangular Top



CENTRAL BASE LOCKING SYSTEM #61777

Easily lock your table into place with our unique system that allows locking from both sides of the table.

Standard Feature: Multi-Specialty

Option: All Ultrasound Tables (Can not be used with Patient Lift Compatible Base)



ULTRASOUND ACCESSORIES

PRONE PILLOW #66180-T01 **For prone or supine positioning!**

The removable pads and low profile design make this the most comfortable positioning solution for your patient. This product offers an added level of comfort for the patient during any ultrasound procedure. Use the pads individually or as part of the system.



Great for Vascular procedures!



T-RAIL ADAPTER #61902

Choose the best spot for this easily installed and removed accessory capable of handling up to 150 lb (68kg). Can be used with any American standard accessory. 8" (20cm) long

Available on all Ultrasound tables.



TABLE EXTENDER #63487-T*

#63487-T*

Lengthens the table by 12" (31cm) to 84" (213cm) allowing for access to the patient's neck for carotid artery exams & thyroid exams. 24" (61cm) W x 12" (31cm) L


Available On All Ultrasound Tables
Requires Accessory Mounting Plate



QUICKLOCK™ HEAD REST & AERO·CEL PAD #PKG6914-T*

The platform allows precise cervical flexion and height positioning over a wide range for both prone and supine treatments.

Available On All Ultrasound tables
Requires Accessory Mounting Plate



DETACHABLE ARM BOARD #62192-T*

Includes clamp to attach Arm Board to the T-rail Adaptor & rotates over 140 degrees. 5" (13cm) wide x 20.5" (52cm) long

Available On All Ultrasound Tables
Requires T-rail Adapter



*Fabric color code must be specified. Fabric options on back cover.

TOP COLORS

TerraTouch Fabric meets the requirements of ISO 10993-1:2009 Biological Evaluation of Medical Tables & passes testing for Cytotoxicity, Skin Irritation & Sensitization as well as CA TB117 Fire Retardance testing. Please refer to the product cleaning instructions for a list of approved cleaners and disinfectants.

T26 White	T09 Butter	T21 Saffron	T12 Opal	T18 Clay	T28 Taupe	T07 Earth
T27 Espresso	T08 Ruby	T19 Orchid	T05 Heron	T20 Ocean	T17 Sapphire	T04 Sky Blue
T30 Seafoam	T22 Sage	T24 Blue Grass	T14 Forest	T13 Stone	T29 Pewter	T01 Coal

We strive to ensure that the colors displayed on our product imagery (online and in print) are as accurate as possible. However, due to the limitations of printers, scanners, computer monitors, and other equipment, the colors may not be accurate.



Oakwork's is independently certificated to ISO 13485 Quality management system and regulatory requirement for medical devices. This registration certifies that Oakworks, Inc.'s quality system conforms to these requirements in our state-of-the-art 90,000 square foot factory in New Freedom, Pennsylvania.



All of our tables and procedure chairs are FDA Listed.

Please visit our website, www.oakworksmed.com, to view our complete line of tables and accessories for C-arm Imaging, Pain Management, Ultrasound, Urology, Lithotripsy, Exam & Treatment.



Call For More Information

Oakworks Medical 717.235.6807 www.oakworksmed.com