

PX SERIES

Exam Table

- 4 designs to choose from
- Up to 3 powered motions
- ADA compliant height range



PX300
shown with Universal Table Extender
and Arm Rest Accessories



PX100



PX150



PX300



ORTHOPEDIC
HILO CASTING
TABLE



ISO 13485:2003
FM 584901

made in the USA
with US & imported parts



OAKWORKS[®]
MEDICAL
Position for Success

PX SERIES

Exam Table

The PX Series Exam Tables give you the flexibility to create the perfect table for your needs at a lower price point. A wide array of options and accessories combine with great ergonomics and patient comfort to deliver a great value.

LOW HEIGHT RANGE FOR EASY PATIENT TRANSFERS

Electric 4
Section Top



MODEL PX300

#70891-T* 31" W 4" Comfort

- Motion 1 Electric Height: 21"-38" (53-97cm)
- Motion 2 Electric Backrest: 0-70°
- Motion 3 Electric Leg Section

Lift-Assist 3
Section Top



ORTHOPEDIC HI-LO CASTING TABLE

#69268-T* 31" W 2.5" Comfort

- Motion 1 Electric Height: 17"-35" (43-89cm)
- Motion 2 Lift-Assist Backrest: 0-90°
- Motion 3 **Manual** Leg Section: 10-90°

Lift-Assist
Backrest Top



MODEL PX150

#70892-T* 31" W 3" Comfort
#72595-T* 31" W 4" Comfort

- Motion 1 Electric Height: 21"-38" (53-97cm)
- Motion 2 Lift-Assist Backrest: 0-80°

Flat Top



MODEL PX100

#70893-T* 31" W 3" Comfort
#73412-T* 31" W 4" Comfort

- Motion 1 Electric 18"-35" (46-89cm) Height

PX SERIES EXAM TABLES

FEATURES	
Top Dimensions - PX100, PX300 Casting	31" (79cm) wide x 73" (185cm) long
Top Dimensions - PX150	31" (79cm) wide x 75" (190cm) long
Foot Control	Controls all powered movements
Weight Capacity	550 lb (250kg)
Product Weight	200 lb (91kg)
Padding - PX100, PX150	3" (8cm) or 4" (10cm) Comfort Foam™
Padding - PX300	4" (10cm) Comfort Foam™
Padding - Casting	2.5" Comfort Foam™
Fabric	Premium, PVC free TerraTouch™
Side Rail	specify side - 23.5" (60cm) long x 6" (15cm) tall (from surface of the table) Orthopedic HiLo Casting Table
Warranty	3 years parts, 2 years labor
Safety Listings	FDA listed, CE marked, cETLus approved, CB Certificate

OPTIONS <i>must be ordered with table.</i>		#
Locking Twin Casters (one end of the table)	To easily move your table in the room	63400
ACCESSORIES <i>order any time.</i>		#
Arm Rests (PX300)	Auto-adjusting & removable for easy access	66760
Universal Table Extender with 3" Padding	21" (53cm) W x 12.5" (32cm) L	68697
Universal Table Extender with 4" Padding	21" (53cm) W x 12.5" (32cm) L	68735
QuickLock™ Head Rest (black) & AeroCel Pad	Comfortable, adjustable prone or supine support	PKG4186-T*
Paper Roll Holder	Vanilla - Accommodates 21" (53cm) roll (roll not included)	11580

*Fabric color code must be specified.

ADJUSTABLE SIDE ARM RESTS #66760-T*

These incredibly versatile, removable arm rests rotate horizontally and always stay level. 7" (18cm) wide at widest point x 23.5" (60cm) long



(PX300 Only)
ACCESSORY

PAPER ROLL HOLDER

#11580 vanilla

This simple design is adjustable to accommodate different size rolls up to 21" (53cm) wide & plugs into the table's outlet holes. *(paper roll not included)*



ACCESSORY

UNIVERSAL TABLE EXTENDER

#68697-T* 3" padding

#68735-T* 4" padding

Perfect as a head rest or foot rest, it has two rods that fit into outlets on the table. Extends table to 84". 21" (53cm) W x 12.5" (32cm) L



ACCESSORY

CASTERS #63400



OPTIONAL

QUICKLOCK™ HEAD REST & AERO·CEL PAD

#PKG4186-T*

The platform allows precise cervical flexion and height positioning over a wide range for both prone and supine treatments.



ACCESSORY

HAND CONTROL

enquire for
part number

Multi-function
hand control



OPTIONAL

TERRATOUGH[™] FABRIC

Meets the requirements of ISO 10993-1:2009 Biological Evaluation of Medical Tables and passes testing for Cytotoxicity, Skin Irritation and Sensitization as well as CA TB117 Fire Retardance testing. For disinfecting purposes you may use Protex, MadaCide, MadaWipe, Accell TB, Virox[®] (Canada), PDI Duper Sani-Cloth Germicidal wipe or a diluted 1/10 bleach solution.



We strive to ensure that the colors displayed on our product imagery (online and in print) are as accurate as possible. However, due to the limitations of printers, scanners, computer monitors, and other equipment, the colors may not be accurate.

PADDING



3" Comfort Foam[™]

*More comfort and support
(Not available on PX 300)*



4" Comfort Foam[™]

Superior comfort and support

CONTACT:



*View our full
line of Medical
products #72703*



Scan this code
with your
smart phone

717.235.6807 www.oakworksmed.com email: customerservice@oakworks.com