

Proning Cushion

Price **\$195**

Part No. 83206



Features:

- Designed for extended proning
- Excellent access to the patient's face
- Accommodates patients from 100lbs to 350lbs
- New improved design has no velcro
- FDA Listed Class I device

Specifications:

DIMENSIONS	27"W x 41"L x 9"H
COVERING	ISO 10993 Compliant
FACE SUPPORT	Adjustable
WEIGHT	7 lbs.
WARRANTY	1 Year



Comfortably prone 100 lbs to 350 lbs patients for extended periods



Recessed area conforms to the patient, cradling and easing patient breathing




Excellent accessibility to the patients face for nasal cannula



OPTIONAL
Storage Case
Part No. 0210
MSRP \$65

made in the USA
with US & imported parts



FM584901

Availability of specific colors is limited.
We reserve the right to substitute other colors to fill your order.

Cover meets the requirements of ISO 10993 Biological Evaluation of Medical devices & passes testing for Cytotoxicity, Skin Irritation & Sensitization.

FDA Listed

Med-MKG-MKT-80 Rev. 2

To place an order, call 717-759-3146 between 8am-5pm EST
or email medsales@oakworks.com