


PX SERIES

The PX line of exam tables feature rock solid support and excellent ergonomics with their twin x beam construction. With a low starting height and a powered range of 18", patient transfers and optimal caregiver ergonomics are facilitated with its barrier free open design. With three top choices and two foam options, the perfect solution for your needs can be found.



PX300
Fully Electric Contour Top

New!
One Piece Removable Pad



Our new design features finger knobs that can be removed in minutes, allowing you to completely remove the entire top and install a new one without any special tools so your table will be in service for much longer.

FEATURES



PX100 Flat Top



PX150 Gas-Assist Backrest



PX300 Fully Electric Contour Top



Integrated Heated Top (Optional)

Our new system is imbedded in the foam, creating a more even heat, and has three settings. Designed originally to be used in car seats and hospital beds, the heating element automatically shuts off after 10 hours providing peace of mind for safety and energy conservation.



SPECIFICATIONS

- New! Removable and replaceable pad
- 3" or 4" Comfort Foam padding
- Large height range - 18" - 35" (PX100)
21" - 38" (PX150 & PX300)

TABLE SPECIFICATIONS

| | |
|--|--|
| Lift System | Electric, Rigid Twin Steel Base Design |
| Controls | Foot Control, Hand Control (Optional) |
| Lifting Capacity | 550 lbs. (250 kg) |
| Fabric | TerraTouch |
| Finish | Powder Coated Steel |
| Padding | 3" (7.6 cm) or 4" (10.2 cm) Comfort Foam (Optional) |
| Motion 1 - Height Range with 3" (7.6 cm) Pad | PX100 - 18" - 35" (46 cm - 91 cm) PX150 & PX300 - 21" x 38" (53 cm - 97 cm) |
| Motion 2 - Backrest | PX300 - Electric Backrest: 0-70° PX150 - Lift Assist Backrest: 0-80° |
| Motion - Leg Rest | PX300 - Electric Leg Rest: 0-30° |
| Width | 29" (73.6 cm), 31" (78.7 cm), 33" (84 cm) |
| Length | 73" (185.4 cm) |
| Options (must be ordered with table) | Integrated Heated Top, Hand Control, Casters |
| Warranty | 3 years parts, 2 years labor |
| Safety Listings | FDA listed, cETLus listed 3034177, CB Certificate |

| | |
|-------------|--|
| Accessories | Adjustable Side Arm Rests, Table Extender, Face Rest Cushion, Quick Lock Face Rest |
|-------------|--|

3" LOCKING CASTERS (OPTIONAL)

In combination with the built in wheels at the other end these casters allow the table to be easily moved in a room.



UNIVERSAL TABLE EXTENDER (ACCESSORY)

Perfect as a head rest or a foot rest, it has two rods that fit into the bracket on the table. Extends table to 85".



ADJUSTABLE SIDE ARMREST (ACCESSORY)

These versatile and comfortable removable arm rests feature a wide range of positions. They always stay level, regardless of the backrest angle. (PX300 only)



TERRATOUGH COLORS

Meets the requirements of ISO 10993-1:2009 Biological Evaluation of Medical Tables & passes testing for Cytotoxicity, Skin Irritation & Sensitization as well as CA TB117 Fire Retardance testing. Please refer to the product cleaning instructions for a list of approved cleaners and disinfectants.



717.759.3146 www.oakworksmed.com email: medsales@oakworks.com